

Disposable Blood Collection Needle



PET Blood Specimen Collection Tube



- ① FPET material has the advantage of safety, low temperature resistance, light weight and environment protection.
- ① Manufactured in the 100,000 grade clean room, free from contamination.
- ① Sterilized by irradiation, which can reduce the use risks.
- ① Accurate blood collection.
- ① High quality rubber, without drop of chips, no precipitate and with good seal property.
- ① There isn't silicon appeared after the barrier gel has been centrifuged. Additive is prayed uniformly on the inner surface of the tube. These properties can protect the precise test machine.
- ① We can professionally custom the blood collection tube for the customer.
- ① We also supply the blood collection needle and needle holder.

Shandong Ande Healthcare Apparatus Co., Ltd.